

VinciSmile®

HugeBond Universal FliPro

8th generation, one component self-etch light cure dental adhesive



➤ Universal Application

- Suitable for various substrates
- Compatible with total-etch, self-etch and selective-etch
- Works with light-cure, self-cure and dual cure cements
- Contains **MDP monomer**, delivering self-etching properties and higher bonding strength

➤ Reliable Bonding

- Bonding strength no lower than 20 MPa
- Low film thickness of 3µm
- Effectively seals over 90% dentin tubules
- Ultra-mild pH(pH=2.7)

➤ Gentle Experience

- Virtually no postoperative sensitivity
- Low occurrence of microleakage prevents secondary caries
- High moisture tolerance offers consistent high bonding strength on dry and moist dentin

*Upgraded flip-open cap
Single-hand operation*

5ML → 6ML



Versatile, Simple, Efficient...




- The latest innovation, **8th generation** adhesive.
- One-bottle technique, streamlines the workflow in practice.
- Featuring universality and versatility, compatible with multiple substrates in restorations.
- Plays a core role in clinical restorative solutions:
 - TopCEM Indirect Restorative Solution
 - TrusFIL Direct Restorative Solution
 - TopCEM Veneer Esthetics Solution




"EASY TO USE WITH MANY INDICATIONS."

Dental Advisor Recommended:

Evaluation Summary:



Compared to Competitive Products:



HUGE

HUGE DENTAL USA LLC
 Web: www.hugedentalusa.com
 Add: 4010 Valley Blvd STE 104, Walnut, CA 91789, USA
 Email: info@hugedentalusa.com Tel: +1 626 283 5808

Follow @HUGE Dental USA